


ACTU8 TALON™

Dynamic Pressure Controlled Forceps




Your Next Generation ACTU8 Forceps!

- Improved visualization
- Less shredding of membranes
- Increased grams of force
- All of the benefits of the ACTU8 handle
- Universal applications including ILM, ERM PVR, dropped IOLs & Diabetic Membranes
- Myopic and standard length available

What is DPC??

Dynamic Pressure Control gives the control of the forceps back to the surgeon. As the surgeon depresses the forceps, it is not an all or nothing with grip strength. Talon tips spread the pressure throughout the body of the tip. This allows for a more gradual build up of grip strength leading to less shredding of membranes and more surgeon control.



TALON TIP PRESSURE AS THE FORCEPS CLOSES

Responsive.

- Linear Actuation
- Instantaneous Tip Movement
- Low Closure Force

Consistent.

- Patented Tip Fabrication Process
- Robust Inspection Procedure
- Reproducible Results, Every Time

Effective.

- Excellent Maneuverability
- Superior Grip Strength¹
- Continuous Evolution





Catalog Numbers and Descriptions

TALON PRO-GRIP²

VS0751.23	23ga Talon Pro-Grip Forceps
VS0751.23-35	23ga Talon Pro-Grip Forceps - Myopic
VS0751.25	25ga Talon Pro-Grip Forceps
VS0751.25-35	25ga Talon Pro-Grip Forceps - Myopic
VS0751.27	27ga Talon Pro-Grip Forceps
VS0751.27-35	27ga Talon Pro-Grip Forceps - Myopic

TALON ILM

VS0752.23	23ga Talon ILM Forceps
VS0752.23-35	23ga Talon ILM Forceps - Myopic
VS0752.25	25ga Talon ILM Forceps
VS0752.25-35	25ga Talon ILM Forceps - Myopic
VS0752.27	27ga Talon ILM Forceps
VS0752.27-35	27ga Talon ILM Forceps - Myopic

TALON END-GRIP

VS0753.23	23ga Talon End-Grip Forceps
VS0753.23-35	23ga Talon End-Grip Forceps - Myopic
VS0753.25	25ga Talon End-Grip Forceps
VS0753.25-35	25ga Talon End-Grip Forceps - Myopic
VS0753.27	27ga Talon End-Grip Forceps
VS0753.27-35	27ga Talon End-Grip Forceps - Myopic

ORDERING INFORMATION:



ACTU8TALON
Dynamic Pressure Controlled Forceps

Vortex Surgical
4 Research Park Drive, Suite 124
Saint Charles, MO 63304 USA

Phone: 636-778-4350
Email: info@vortextsurgical.com
Web: www.vortextsurgical.com

1-Vortex Data on file 2-Coming Soon



VORTEX SURGICAL