

SURGICAL LENSES

Available in 6 popular styles, these lenses deliver high resolution direct-image retinal views for all vitrectomy procedures. Most are fitted with a silicone stabilizing ring, eliminating the need for a suture ring or other lens holding device.



BICONCAVE

PRODUCT# VBCD

- Direct image for fluid air exchange
- 25° field of view
- 0.8x magnification
- Sterile, 10/Box



FLAT

PRODUCT# VFD

- Routine direct image of central retina
- 36° field of view
- 1.0x magnification
- Sterile, 10/Box



FLAT SSV

PRODUCT# VFLATSSVD

- Routine direct image of central retina
- 36° field of view
- 0.92x magnification
- Sterile, 10/Box



MAGNIFYING

PRODUCT# VMD

- High magnification of central retina
- 30° field of view
- 1.5x magnification
- Sterile, 10/Box



30° PRISM

PRODUCT# V30PD

- Off axis direct image
- 33° field of view
- 1.0x magnification
- Sterile, 10/Box



WIDEFIELD

PRODUCT# VWFD

- Wide field direct image
- 48° field of view
- 0.5x magnification
- Sterile, 10/Box

