

I² KITS

FOR INTRAVITREAL INJECTIONS

Injection kits are comprised of the tools you need to perform sterile intravitreal injections in the office or operating room.

We offer three popular configurations, however, we would be happy to design a kit to meet your needs.



BASIC

- Adult Blade Speculum
- Scleral Marker



STANDARD

- Adult Blade Speculum
- Scleral Marker
- 1cc Syringe (2)
- 18ga x 1" Hypo Needle (2)
- 27ga x 1/2" Hypo Needle
- 30ga x 1/2" Hypo Needle
- Cotton Tip Applicator (4)
- Exam Gloves
- Soaking Tray
- Small Mayo Tray Cover

STANDARD+DRAPE

- Adult Blade Speculum
- Scleral Marker
- 1cc Syringe (2)
- 18ga x 1" Hypo Needle (2)
- 27ga x 1/2" Hypo Needle
- 30ga x 1/2" Hypo Needle
- Cotton Tip Applicator (4)
- Exam Gloves
- Soaking Tray
- Small Mayo Tray Cover
- Fenestrated Eye Drape

Catalog Numbering and Ordering information are located on the reverse side





I2 KITS FOR INTRAVITREAL INJECTION

VS0500	I2 Injection Kit - Basic, Sterile, 10/bx
VS0510	I2 Injection Kit - Standard, Sterile, 10/bx
VS0520	I2 Injection Kit - Standard + Drape, Sterile, 10/bx

**To design your own custom I2 Injection Kit,
call - 636-778-4350 or
email - info@vortexsurgical.com**

ORDERING INFORMATION:



I2 INJECTION KITS

Vortex Surgical
4 Research Park Drive, Suite 124
Saint Charles, MO 63304 USA

Phone: 636-778-4350
Email: info@vortexsurgical.com
Web: www.vortexsurgical.com



VORTEX SURGICAL