




YOUR PARTNER IN OCULAR ONCOLOGY

Vortex Surgical prides itself on listening and supporting the needs of today's Ocular Oncologists. From new illumination devices to ground breaking biopsy devices, Vortex Surgical is here for you.



BEACON ILLUMINATION SYSTEM



BEACON LED LIGHT SOURCE

- 50-watt LED Light Source
- Inexpensive & Portable
- Perfect for the Clinic or the Operating Room



BEACON TRANSILLUMINATOR

- One-Piece Transilluminator
- Side-Fire Illumination allows greater access
- Clinic or Operating room use



TRANSPUPILLARY ILLUMINATOR

- Atraumatic Cornea Shield & Bright, Controlled Illumination
- Tumor identification & Marking
- Clinic or Operating Room Use

NANO SUBRET GATEWAY DEVICE



- 28ga Hypodermic needle with 41ga polyimide biopsy needle
- Perfect for tumor biopsy
- Curved and Straight versions available
- No need for vitrectomy
- Clinic or Operating Room Use

EDD - EXTERNAL DRAINAGE & DEPRESSION



- 28ga hypodermic needle with 2.4 mm extension
- Works well for different types of biopsies
- Curved and Straight versions available
- No need for vitrectomy
- Clinic or Operating Room Use



Catalog Numbers and Descriptions

BEACON ILLUMINATION SYSTEM

VS0900	Beacon LED Light Source, 1/bx
VS0801	Transilluminator for Alcon, B&L and DORC Vitrectomy systems, Sterile, 5/bx (Adapter needed)
VS0802	Beacon Transilluminator for Beacon LED Light Source Sterile, 5/bx
VS0803	Transpupillary Illuminator for Beacon Illumination System Sterile, 5/bx
VS0804	Transpupillary Illuminator for Alcon, B&L and DORC Vitrectomy systems, Sterile, 5/bx (Adapter needed)

BIOPSY OPTIONS

VS0225	Straight NANO SubRet Gateway Device, Sterile, 5/bx
VS0225C	Curved NANO SubRet Gateway Device, Sterile, 5/bx
VS0290	Straight EDD - External Drainage & Depression Device, Sterile, 5/bx
VS0290C	Curved EDD - External Drainage & Depression Device, Sterile, 5/bx

ILLUMINATION ADAPTERS

VS0ALC	Illumination Adapter - Alcon Constellation and Alcon Accurus, 1/bx
VS0BLS	Illumination Adapter - B&L Stellaris and B&L Stellaris PC, 1/bx
VS0DRC	Illumination Adapter - DORC EVA and DORC NEXUS, 1/bx

ORDERING INFORMATION:



Vortex Surgical Website

Vortex Surgical
4 Research Park Drive, Suite 124
Saint Charles, MO 63304 USA

Phone: 636-778-4350
Email: info@vortexsurgical.com
Web: www.vortexsurgical.com



Vortex Ocular Oncology